Guia do utilizador do Instron Connect - Menu principal

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A screenshot of the Instron Connect sign in/sign up form. The form is titled 'INSTRON The difference is measurable'. It contains a sign in section with an 'Email Address' input field and a 'Password' input field, with a 'Forgot your password?' link. Below this is a blue 'Sign in' button. At the bottom, there is a link for 'Don't have an account? Sign up now' where 'Sign up now' is highlighted with a red box.

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INSTRON
The difference is measurable

Email Address
Email Address
Send verification code

New Password
New Password

Confirm New Password
Confirm New Password

Given Name
Given Name

Surname
Surname

Display Name
Display Name

Create Cancel

INSTRON Instron Connect HOME AGREEMENT SOFTWARE SUPPORT CERTIFICATES SERVICE HISTORY TROUBLESHOOTING ADMIN English BELLA

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
The screenshot displays the 'Support Request' form. At the top, the 'Request Type' is set to 'Instron Connect Portal Access'. Below this, the 'REQUESTOR' section contains fields for 'First Name', 'Last Name', 'Email Address', and 'Phone number'. The 'Phone number' field is highlighted with a red box. Underneath, there are radio buttons for 'Preferred Contact', with 'Email' selected. The 'ADDITIONAL USER' section follows, with fields for 'First Name', 'Last Name', 'Email Address', and 'Phone number'. Below these fields is a dropdown menu for 'Agreement Number' set to '--Contract Number--'. A large text area for 'Message' is provided at the bottom, and a red 'Submit' button is located at the very bottom of the form.

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My Agreement

Agreement Information:

Agreement Number	From Date	End Date
RTS: [redacted]	01/24/2022	01/23/2023

Equipment covered under the selected agreement:

Serial Number	System ID	Customer Asset Number	Coverage Type
3345 [redacted]	3345 [redacted]		Instron Connect
34TM [redacted]	34TM [redacted]		Instron Connect
68TM5C [redacted]	68TM [redacted]		Instron Connect

RTS0022 includes the following scope of work:

Serial Number	Description	Service Description	Qty
Machine: 3345 [redacted] - Tag #: 3345 [redacted]			
3345 [redacted]	5 kN Single Column	RTS Agreement powered by Connect (Multiple Systems)	1
Machine: 34TM [redacted] - #: 34TM [redacted]			
34TM [redacted]	5KN, 3400 SERIES, TABLE MODEL LOAD FRAME	RTS Agreement powered by Connect (Multiple Systems)	1
Machine: 68TM50 [redacted] - Tag #: 68TM50 [redacted]			
68TM50B [redacted]	50KN, 6800 SERIES, TABLE MODEL LOAD FRAME	RTS Agreement powered by Connect (Multiple Systems)	1

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Instron Connect Agreements include the following features and benefits:

6

Instron Technical Support Services:

- Calibration and Contract Renewal Reminders
- Software Update Notifications
- Direct Downloadable Software Updates
- Downloadable Software Updates via Instron Connect Portal
- Calibration Certificates, Service History and Contract Information via Instron Connect Portal
- System Diagnostics via Instron Connect Software
- Direct Screen Share with Voice through Instron Connect Software
- Augmented Reality (AR) App Support
- Self-Diagnostics Tool via InSkill Mobile App and Instron Connect Portal

RTS includes the following scope of work:

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
Software Updates

Use this page to download a selected software.

*Please note: If your system is IQ/OQ validated for compliance with FDA (or equivalent) standards, do not install or update Instron software on your PC without prior approval from your internal Quality Team.

If you believe you should have access to download additional software please click [here](#)

System ID: Please select a Software Type:

File Name	Download File
IC_2.02.8100.zip	

Software Updates

Use this page to download a selected software.

*Please note: If your system is IQ/OQ validated for compliance with FDA (or equivalent) standards, do not install or update Instron software on your PC without prior approval from your internal Quality Team.

If you believe you should have access to download additional software please click [here](#)

System ID: Please select a Software Type:

File Name	Download File
BHU_4.34.30245.zip	

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
34TM

Message

Attach Files Choose Files No file chosen

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INSTRON Instron Connect

HOME AGREEMENT SOFTWARE SUPPORT CERTIFICATES SERVICE HISTORY INSKILL

Successful Request Submission

Your request has been submitted successfully. Please verify you have received an email regarding your submission.

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Certificates

To view PDF file, you must have [Acrobat Reader](#) installed on your system.
If your certificates do not appear, please [submit a helpdesk ticket](#)

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CERTIFICATE OF CALIBRATION
ISSUED BY: INSTRON CALIBRATION LABORATORY
DATE OF ISSUE: 26-Aug-2022
CERTIFICATE NUMBER: 21208232092909

INSTRON
825 University Avenue
Norwood, MA 02062-2643
Telephone: (800) 473-7838
Fax: (781) 875-9750
Email: service_requests@instron.com

NVLAP
CALIBRATION
www.nvlap.com/2020/04/

Page 1 of 5 pages
APPROVED SIGNATORY

Verification Results:
System ID: 34TM
Indicator - Service Port (in)
Range: 0.10043 in to 4.00581 in - Ascending Starting Position: 20.00 in Maximum Error: -0.213%
Range: -0.10024 in to -4.01083 in - Descending Starting Position: 24.00 in Maximum Error: -0.746%
PASSED Class A
The starting position is measured from the base seen in the bottom of the crosshead.

Customer:
Name: Instron Training Center
Address: 825 University Avenue, Norwood, MA 02062, United States
Contact: Joseph Curaso, joseph_curaso@instron.com
Service Order No.:

Machine/System:
Manufacturer: Instron
Serial No.: 34TM
Condition: Good

Temperature:
Starting Temperature: 71.4 °F
Final Temperature: 71.4 °F

Methodology:
The assessment of the testing machine was conducted on site at the above customer location in accordance with ASTM E2309/E2309M-20 "Standard Practice for Verification of Displacement Measuring Systems and Devices Used in Material Testing Machines" (Follow-the-Displacement Method) using Instron procedure KA-8-07.
The system was calibrated in the 'As Found' condition with no adjustments or repairs carried out. This is also the 'As Left' condition.
Prior to verification, a pre-calibration inspection was conducted. During the inspection, the testing system was found to be in Good condition.
The calibrated range of displacement includes only those displacements which are greater than or equal to the ASTM Lower Limit.
These calibration tests were made with the testing machine in the vertical position.

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Service History

To view PDF file, you must have [Acrobat Reader](#) installed on your system.

If your service history does not appear, please [submit a helpdesk ticket](#)

System ID: *

- Select System ID--
- 3345
- 34TM
- 68TM50

Order ID	Status	Service Type	Service Description	Created Date	Last Updated Date	View Service Report
<input type="text" value="Filter..."/>	<input type="text" value="Filter..."/>	<input type="text" value="Filter..."/>	<input type="text" value="Filter..."/>	<input type="text" value="Filter..."/>	<input type="text" value="Filter..."/>	
SV21	Closed	PM/Calibration	training calibration call	May 11, 2021	May 24, 2021	
SV21	Closed	PM/Calibration	calibration training call	May 11, 2021	May 19, 2021	

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System ID: 34TM

Order ID	Status	Service Type	Service Description	Created Date	Last Updated Date	View Service Report
Filter...	Filter...	Filter...	Filter...	Filter...	Filter...	
SV21	Closed	PM/Calibration	training calibration call	May 11, 2021	May 24, 2021	
SV21	Closed	PM/Calibration	calibration training call	May 11, 2021	May 19, 2021	

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Company: Instron Training Center
Address: 825 University Avenue
Norwood MA 02062
United States

Contact:
Service Order #: SV21
Service Engineer:
Purchase Order #:

Machine Serial #: 34TM-
Machine Tag #: 34TM-
Customer Asset #:

Labor Activities*

Date	Qty	Product #	Product Description
5/11/2021	1	S1460-201-A	ENH Speed & Displacement Verification (Both directions)
5/12/2021	1	S1481-202-A	ENH Force CAL 500N-100KN (both directions)
5/21/2021	1	S1490-802-A-ES-A	ASTM Strain ENH Single Dir for Clip-On, XL,COD,LVDT


*Note: Labor activities only reflect the type of work performed on this date and not the actual hours worked.

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

A screenshot of a web interface for a diagnostic tool. It features three dropdown menus. The first two are set to '34TM' and '2712-04x Pneumatic Grips'. The third dropdown menu is open, showing 'Diagnose' as the selected option, with 'Diagnose' and 'How Do I?' as visible options in the list.

A screenshot of the same web interface. The dropdown menus are now closed. A red 'Begin' button is highlighted with a red border, indicating it is the next step in the process.

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
A screenshot of a web interface titled 'What error are you experiencing?'. It features a dropdown menu with 'SELECT ONE' at the top. The menu is open, showing a list of error types: 'Run time diagnostic failure' (highlighted in blue), 'Self test diagnostic failure', 'Air kit problem', 'Load cell won't calibrate', 'Communications error', and 'Bluehill Software Issue'.

A screenshot of a web interface titled 'What symptoms are you experiencing?'. It features a dropdown menu with 'Grip closing slowly' selected. Below the dropdown are two buttons: a blue 'Submit' button and a white 'Skip' button with a right-pointing arrow.

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


Please select your country or region

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Take a photo that helps illustrate the issue

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Please provide any additional information that may help diagnose your problem.

Grips are still closing slowly, we have tried adjusting

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 **Your support case has been sent**

Your request has been forwarded and a representative will respond to you as soon as possible.
If you have any further questions, you can reach us at 800-473-7838, press 3 for service support.

1 

[Find your local office on our locations page to call for tech support.](#)

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
Support Request

Calibration

Submit Request
Enter your information below then click submit

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


ID	Time	Message
600	30/06/2022 12:21:10	Remote access is now enabled.
105	30/06/2022 12:19:15	Successfully connected to Instron® Connect.


Partilhar o ecrã através do software Instron Connect (2/3)

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

ID	Time	Message
605	30/06/2022 12:24:48	Remote access is now disabled.
600	30/06/2022 12:21:10	Remote access is now enabled.
105	30/06/2022 12:19:15	Successfully connected to Instron® Connect.

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
Message Center		
ID	Time	Message
105	30/06/2022 14:36:32	Successfully connected to Instron® Connect.
315	30/06/2022 14:36:32	Calibration was performed over 12 months ago. If you wish to recalibrate your system, please contact Instron Service under the 'Calibration' tab to schedule a visit.

System Status		
	Cloud Connection	
	Calibration	
<input type="button" value="Enable"/>	Remote Access	

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To activate your account, click on the button below, and use the activation code when prompted to start using InSkill. To get started, check out [InSkill University](#), found under your user icon in the **InSkill Web Application**.

1 **Code:** [blurred]

2 **Setup Account**

[Click here to setup your account](#)

Smart Products Deserve Smart Support

Welcome to InSkill

Please enter the code in your invite email to verify your identity and activate your account.

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
A screenshot of the 'Team' page in the application. The 'Team' menu item is highlighted with a red box. Below the navigation bar, there is an 'Export' button and a filter dropdown set to 'Active'. The main content is a table with columns for Name, Role, and Copilots.

Name	Role	Copilots
bob_jones	Helpdesk Agent	All
christine_smith	User 1 - InSkill Apps	All

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

Help Escalate ×

Request specific content types:

- *Show me videos about maintenance*
- *Find troubleshooting tasks*

Commands

- **Save solution:** ... Saves a new solution
- **Save note:** ... Saves a new note to the asset (when asset is selected)

Troubleshooting

"I could not find an answer":

- The information isn't in your indexed documentation
- Check document availability: Ensure relevant documents are indexed

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Step History Data Resources Handoff Escalate powered by InSkill

We were unable to find a solution for your issue. The next step is to escalate your diagnosis to InStron Support. Please select the button below to start that process.


Continue...

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
Passo seguinte


1. Selecione Gerir > InSkill


2. Selecionar tarefas.

3. Select General Utilities as a Copilot. Then, Select Contact Support as a Task.

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First Name: Harold

Preferred Contact: Email Phone

- Select Request Type--
- Instron Connect Portal Access
- Technical Support
- Calibration Quote
- On-site Service
- Service Parts or Returns
- Other

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